SPHEROTM Calibration Particles

SPHEROTM Calibration Particles are designed for routine calibration of flow cytometers. They are used extensively by many laboratories for QC and long term performance tracking of the flow cytometer. In addition, they are also used for routine alignment and calibration in fluorescence and confocal fluorescence microscopy.

SPHERO[™] Rainbow Calibration Particles

The **Rainbow Calibration Particles** (RCPs) contain a mixture of several similar size particles with different fluorescence intensities. Every particle contains a mixture of fluorophores that allows excitation at any wavelength from 365 to 650 nm. As a result, all channels in the flow cytometer can be calibrated using the same set of particles.

The RCPs provide a reliable and reproducible particle mixture for calibrating flow cytometers. They are very stable since the fluorochromes are entrapped within the particles instead of being located on the surface. In addition, Spherotech uses fluorophores that are non-spectral matching to the commonly used fluorophores such as FITC, PE or PE-Cy5. As a result, the RCPs are stable in terms of fluorescence.

The RCPs are convenient and affordable to use for long term performance tracking. They are packaged in a dropper bottle to facilitate the dispensing and storage. The diluted particles can be stored in the freezer for later use if desired to reduce costs. Dilution of a few drops of the particles from the dropper bottle to 1 mL of a diluent will provide adequate particle concentration for flow cytometer calibration. The diluted Rainbow Calibration Particles remain stable following repeated freezing and thawing.

Particle Type and Surface	Size, µm	Catalog No.	Unit
Rainbow Calibration, 4 peaks, 107/mL	1.8-2.2	RCP-20-5	5 mL
Rainbow Calibration, 6 peaks, 107/mL	3.0-3.4	RCP-30-5	5 mL
Rainbow Calibration, 8 peaks, 107/mL	3.0-3.4	RCP-30-5A	5 mL
Rainbow Calibration, 4 peaks, 107/mL	3.5-4.0	RCP-35-5	5 mL
Rainbow Calibration, 6 peaks, 107/mL	6.0-6.4	RCP-60-5	5 mL

The contour plot below displays the bead distributions of the Rainbow Calibration Particles in FITC vs PE. (Catalog No. RCP-30-5A).



Due to the stability of **SPHERO**[™] Rainbow Calibration Particles, these particles are very useful for long term tracking of instrument performance as well as for routine calibration of the flow cytometer. Advantages of using water insoluble fluorescent dyes to obtain fluorescent particles as compared to surface labeled microparticles are listed below.

The **SPHERO**[™] Fluorescent Particles and Calibration Particles:

- Stable for several years when stored properly.
- Withstand freeze-thaw cycles; diluted particles can be stored frozen for later use.
- Contains multiple fluorophores incorporated in the same particle to be used in different channels of the flow cytometer.
- Available with different fluorescent intensities on the same size particles.
- Can be sanitized by treating with 70% ethanol or other antibiotic agents.

Histograms showing individual peaks representing various fluorescence intensities in Rainbow Calibration Particles (Catalog # RCP-30-5A) are shown below.



Tel.: 800-368 0822 or 847-680 8922; Fax: 847-680 8927; E-Mail: service@spherotech.com Visit us on the web at *http://www.spherotech.com*









Tel.: 800-368 0822 or 847-680 8922; Fax: 847-680 8927; E-Mail: service@spherotech.com Visit us on the web at *http://www.spherotech.com* The relative number of fluorophores per particles has been determined for every peak of RCP-30-5 in FL1 (FITC, MEFL), FL2 (RPE, MEPE), FL3 (RPE-Cy5, MEPCY) and FL4 (APC, MEAP) channels of flow cytometer to plot the calibration graph as shown below. The calibration graph is used to check the linearity of the PMT in each channel. In addition, the relative number of fluorophores can be cross calibrated with cells or particles stained with known number of spectral matching fluorophores such as FITC, PE, RPE-Cy5 to estimate the number of fluorophores on stained cells. The RCP-30-5A, which is identical to RCP-30-5 with the exception of two additional peaks between the blank and the dimmest peak of RCP-30-5 to give a total of 8 peaks is shown on Page 16. The RCP-30-5A is very useful in checking the sensitivity and resolution of the flow cytometer.



A Template for MS Excel files, as shown below, is free of charge upon request. The template will allow the user to check and report the linearity of PMT in all channels easily by using RCP-30-5, RCP-30-5A, RCP-60-5, URCP-38-2K, URCP-50-2K, RQC-4K or ACP-30-5K.



The Rainbow Calibration Particles are available in sizes of 2 to 6 µm as shown below to suit various applications.



Tel.: 800-368 0822 or 847-680 8922; Fax: 847-680 8927; E-Mail: service@spherotech.com Visit us on the web at *http://www.spherotech.com*