

## **Magnetic Particle Coated with Pepsin and Papain**

### **Preparation of IgG F(ab')<sub>2</sub> Fraction Using Pepsin Coated Magnetic Particles (Spherotech Cat. # PEPM-40-2)**

The Spherotech Cat. No. PEPM-40 consists of 4µm magnetic particles covalently linked to pepsin from porcine stomach mucosa. These particles efficiently cleave IgG from various species producing F(ab')<sub>2</sub> and other lower molecular weight products within a few hours. After digestion Pepsin-Magnetic particles can be easily removed magnetically from the reaction vessel leaving no pepsin in the solution. As a result, the cleavage of IgG can be effectively controlled. The supernatant containing products from the cleavage and any residual IgG can be further purified chromatographically or by using Protein A magnetic particles (Spherotech Cat. # PAMX-10).

This procedure allows the controlled preparation of essentially pure F(ab')<sub>2</sub> in simple two steps. In addition, pepsin magnetic particles can be reused repeatedly without a significant loss in activity.

#### **Protocol:**

##### **Materials Needed:**

1. 0.01M acetate buffer, pH 3.0
2. Spherotech Flexi-Mag Magnetic Separator (Spherotech Cat. # FMJ-1000)
3. Incubator set at 37 °C.
4. Rotator or stirrer.
5. Protein A coated magnetic particles (Spherotech Cat # PAMX-10)

#### **Procedure:**

1. Make an IgG solution in 0.01M acetate buffer, pH 3.0. The preferred concentration is 1mg/mL.
2. Transfer 200-500µL of Pepsin-magnetic particles (1% w/v) for each mg of IgG to a screw capped vial or tube.
3. Separate the particles magnetically. Remove the supernatant. Wash the particles 1 time with acetate buffer, pH 3.0.
4. Add IgG solution in acetate buffer, pH 3.0 to the washed Pepsin-magnetic particle pellet.
5. Cap, mix and rotate or stir gently for 2 to 24 hours at room temperature or 37°C. The duration of incubation can be adjusted to achieve the desired level of IgG cleavage.
6. Magnetically separate the particles and collect the supernatant. Transfer the particles to a vessel containing the neutralizing solution ( e.g. ~ 30µL of 10N NaOH for each mL of solution).
7. Add Protein A magnetic particles (PAMX-10) to the supernatant from step 6 to bind and remove residual IgG and Fc fractions from the solution.
8. The supernatant remaining after PAMX-10 treatment will contain mostly F(ab')<sub>2</sub> fragments.

#### **Note:**

- Each mL of PAMX-10 will bind to 30-50µg of human, mouse and rabbit IgG; however, IgG of some species such as rat and goat will not bind to protein A.
- Collect the Pepsin magnetic particles and Protein A magnetic particles (PAMX-10) for reuse. Wash particle pellet with PBS or other neutral buffer and store refrigerated.

### **Preparation of IgG Fab Fraction Using Papain Coated Magnetic Particles (Spherotech Cat. # PAPM-40-2)**

The Spherotech PAPM-40 consists of magnetic particles (4 $\mu$ m) covalently linked to Papain from Papaya latex. These particles efficiently cleave IgG into Fab and Fc fragments within a few hours. After digestion Papain-Magnetic particles can be easily removed from the reaction vessel magnetically leaving no Papain in the solution. As a result, effectively stopping further cleavage of IgG. The supernatant containing products of cleavage and any undigested IgG can be further purified chromatographically.

This procedure allows the controlled preparation of essentially pure Fab in simple two steps. In addition, papain magnetic particles can be reused repeatedly without a significant loss in activity.

#### **Protocol:**

Materials Needed:

1. 0.01M acetate buffer, pH 5.0
2. 20mM NaEDTA solution in water
3. 1M L-Cysteine solution in 1N HCl
4. 0.5M Iodoacetamide solution in water
5. Spherotech Flexi-Mag Magnetic Separator (Spherotech Cat. # FMJ-1000)
6. Incubator set at 37°C.
7. Rotator or stirrer.

#### **Procedure:**

1. Make an IgG solution in 0.01M acetate buffer, pH 5.0. Add 50 $\mu$ L of EDTA and 25 $\mu$ L of L-Cysteine solutions to each mL of IgG solution. Preferred concentration of IgG is 1mg/mL.
2. Transfer 200-500 $\mu$ L of Papain-magnetic particles (1% w/v) for each mg of IgG to a screw capped vial or tube.
3. Separate the particles magnetically. Remove the supernatant. Wash the particles 1 time with acetate buffer, pH 5.0.
4. Add the IgG solution in acetate buffer, pH 5.0 to the washed Pepsin-magnetic particle pellet.
5. Cap and mix by rotation or stirring gently for 2 to 24 hours at room temperature or 37°C. The duration of incubation can be adjusted to achieve the desired level of IgG cleavage.
6. Magnetically separate the particles and collect the supernatant. Add 100  $\mu$ L of Iodoacetamide solution to each mL of supernatant to stabilize the Fab.

**Note:** Collect the Papain-magnetic particles for reuse. Wash the particle pellet with PBS or other neutral buffer. Refrigerated the particles during storage.

#### **Why Magnetic Particle Coated Enzymes ?**

Several advantages of using enzymes covalently linked to magnetic particles over enzymes attached to non magnetic gel, latex, or soluble enzymes are listed below:

1. The magnetic separation of enzymes covalently linked to magnetic particles stops the cleavage of IgG without additional inhibitors, pH changes or other extraneous agents. This avoids contamination of the sample.
2. Magnetic particles coated with enzymes are can be reused repeatedly without any significant loss in enzymatic activity.
3. The reaction kinetics of small uniform magnetic particles coated with various enzymes is fast and efficient.

---

**BIBLIOGRAPHY-Selected references**

1. Arend WP, et al. Characterization of two species of Fc fragments obtained from rat IgG2a by prolonged papain digestion. *J Immunol.* 1977 Feb; 118(2): 388-394.
2. Butler JE, et al. The differential enzyme susceptibility of bovine immunoglobulin G1 and immunoglobulin G2 to pepsin and papain. *Biochim Biophys Acta.* 1978 Jul 21; 535(1): 125-137.
3. Harlow Ed and Lane D. In *Antibodies-A laboratory Manual.* By Cold Spring Harbor Laboratory. 1988
4. Kojima M, et al. Amino acid sequence of the variable region of heavy chain in immunoglobulin (Mot) having unusual papain cleavage sites. *Mol Immunol.* 1986 Feb; 23(2): 169-174.
5. Kojima M, et al. An unusual papain cleavage of a human IgG1 (lambda) myeloma protein (Mot). *Mol Immunol.* 1982 Sep; 19(9): 1095-1103.
6. Nardella FA, et al. Fc intermediate (Fci), a papain-generated fragment of human IgG, intermediate in charge, molecular weight and cleavage between the Fc and Fc' fragments of IgG. *Mol Immunol.* 1985 Jun; 22(6): 705-713.
7. Porter R R, et al. Hydrolysis of rabbit gamma globulin and antibodies with crystallin papain. *Biochem J* 1959; 73: 119.
8. Rousseaux J, et al. Optimal conditions for the preparation of Fab and F(ab')<sub>2</sub> fragments from monoclonal IgG of different rat IgG subclasses. *J Immunol Methods.* 1983 Nov 11; 64(1-2): 141-146.
9. Stewart GA, et al. The effect of acid treatment upon the susceptibility of rabbit IgG to proteolytic cleavage with various enzymes. *Immunochemistry.* 1975 Aug; 12(8): 713-721.
10. Steinbuch M, et al. Proteolytic digestion and chemical cleavage of human IgA. *Prog Immunobiol Stand.* 1970; 4: 77-80.
11. Tseng J, et al. Differential susceptibility of rabbit IgA subclasses to trypsin cleavage: characterization of the fragments obtained from the g subclass. *J Immunol.* 1975 Aug; 115(2): 454-461.
12. Utsumi S. Stepwise cleavage of rabbit immunoglobulin G by papain and isolation of four types of biologically active Fc fragments. *Biochem J.* 1969 Apr; 112(3): 343-355.
13. Wang AC, et al. Cleavage sites of human IgG1 immunoglobulin by papain. *Immunochemistry.* 1977 Mar; 14(3): 197-200.
14. Zappacosta S, et al. Mechanism of cleavage of rabbit IgG in two stages by soluble papain and reducing agent. *J Immunol.* 1968 Jun; 100(6): 1268-1276.