


**SPHERO™ Ultra Rainbow Calibration Particle Kits**

A new set of calibration particles, **UltraRainbow Calibration Particles Kits**, are now available for performance tracking of flow cytometers with fluorescent channels in the Far Red. The **UltraRainbow Calibration Particles** are available in either 3.8µm or 5.1µm. The **UltraRainbow Calibration Particles** have improved resolution in the PE-Cy7, APC, and APC-Cy7 channels. These particles are supplied as a kit. One 3 mL bottle of blank particles and a 2 mL bottle of fluorescent particles.

Particle Type and Surface	Size, µm	Catalog No.	Unit
UltraRainbow Calibration, 6 peaks, 10 <sup>7</sup> /mL	3.5-3.9	URCP-38-2K	1x3mL & 1x2 mL
UltraRainbow Calibration, 6 peaks, 10 <sup>7</sup> /mL	5.0-5.9	URCP-50-2K	1x3mL & 1x2 mL

**Histograms of URCP-50-2K**



**Histograms of URCP-38-2K**

