

Spherotech, Inc.


27845 Irma Lee Circle, Unit 101, Lake Forest, IL 60045

SPHERO™ Ultra Rainbow Calibration Particle Kit


A new set of calibration particles, **Ultra Rainbow Calibration Particles**, are now available for performance tracking of flow cytometers with fluorescent channels in the Far Red. The **Ultra Rainbow Calibration Particle Kits** are available in either 6 peaks, 3.8µm or 5.05µm. These particles are supplied as a kit. One 3mL bottle of blank particles and a 2 mL bottle of fluorescent particles. The **Ultra Rainbow Calibration Particles** have improved resolution in the PE-Cy7, APC, and APC-Cy7 channels.

Particle Type and Surface	Size, µm	Catalog No.	Unit
Ultra Rainbow Calibration Kit, 6 peaks, 10 ⁷ /mL	3.5-3.9	URCP-38-2K	2 mL
Ultra Rainbow Calibration Kit, 6 peaks, 10 ⁷ /mL	5.0-5.4	URCP-50-2K	2 mL

Histograms of URCP-38-2K



Histograms of URCP-50-2K



Tel.: 800-368 0822 or 847-680 8922; Fax: 847-680 8927; E-Mail: service@spherotech.com

Visit us on the web at <http://www.spherotech.com>