


SPHERO™ Ultra Rainbow Fluorescent Particles for Beckman Coulter Cytometers

- Consists of a single peak for performance verification of Beckman Coulter flow cytometers
- Determines cleanliness of the flow cell and fluidics

Particle Type and Surface	Size, μm	Catalog No.	Unit
Calibration Beads for MoFlo™ Astrios™, $10^7/\text{mL}$	3.0-3.4	B27034	2 mL
Calibration Beads for MoFlo™, MoFlo™ XDP, and CyAn™ ADP, $10^7/\text{mL}$	3.0-3.4	B28479	2 mL
AccuCount Ultra Rainbow Fluorescent Particles (Calibration Beads for Aquios CL™, $10^6/\text{mL}$)	3.8 (+/-0.3)	ACURFP-38-5	5 mL
AccuCount Ultra Rainbow Fluorescent Particles (Calibration Beads for Aquios CL™, $10^6/\text{mL}$)	3.8 (+/-0.3)	ACURFP-38-15	15 mL



Dot Plot and histograms of B28479 on a Beckman Coulter Cyan™ ADP